
## QuikPath's SARS CoV-2 Isothermal RT-PCR





## Easy-to-Use by Anyone, Anywhere!

Insert oral / nasal swab into cartridge to activate



Place the catrdige in heat block



Visual result ready in 40 mins



Scan barcode on cartridge using QuikPath app to store results



Imal RT. PCP

## Advanced Technology, Easy Storage


- Based on LAMP RT-PCR technology, our tests can detect viral loads as low as 10 viral copies/µL
- Microfluidic cartridge design eliminates need for pipettes and skilled personnel
- Contains internal control
- Multiplexed 2-gene assay for detection of multiple strains
- Our test can detect new UK and South African COVID-19











220 Westwood Avenue, #11-06, The Floravale, SG 648352

For pricing and further info, contact us:

QuikPath Pte. Ltd. (Reg. No. 201937852G)







## info@quikpath.sg

\* As determined through lab-based analytical studies \*\* As determined via in-silico testing of NCIB and GISAID databases